



System Requirements Healthquest 4.0

Minimum Requirements for workstations:

- Dual-Core, i-Series or Better
- 500GB hard drive
- 19-inch monitor, at 1024 x 768 resolution
- Windows 7 Pro or better (32 bit or 64 bit), or Windows 8/8.1
- 4GB of RAM
- High speed Internet account (ADSL or Cable if available)
- Twain drivers for scanning
- Microsoft .NET Framework 3.5 SP1

Recommended Requirements for workstations:

- Quad-Core, i3 or Better
- 128GB SSD hard drive
- 24-inch monitor, at 1080p resolution
- Windows 7 Pro or better (64 bit), or Windows 8/8.1
- 8GB of RAM
- High speed Internet account (ADSL or Cable if available)
- Twain drivers for scanning
- Windows compliant label printer (Dymo recommended)
- Laser Printer
- Microsoft .NET Framework 3.5 SP1

Minimum Requirements for servers:

- Dual-Core, i-Series or Better
- 1TB HDD RAID compliant
- 19-inch monitor, at 1024 x 768 resolution
- Windows Server 2008 or better
- 4 GB ECC RAM
- High speed Internet account (ADSL or Cable if available)
- Twain drivers for scanning
- External drive or some type of backup drive
- Microsoft .NET Framework 3.5 SP1

Recommended Requirements for servers:

- Quad-Core, i-Series or Better (64-bit)
- 1TB+ HDD RAID hybrid, or SSD
- 24-inch monitor, at 1080p resolution
- Windows Server 2008 or better
- 8 GB ECC RAM
- High speed Internet account (ADSL or Cable if available)
- Twain drivers for scanning
- External drive or some type of backup drive
- Microsoft .NET Framework 3.5 SP1

Setup Requirements Prior to Installation

- Electronic Faxing Software for incoming and outgoing faxes.
- All users setup with windows password and permissions.
- Backup plan in place.
- Windows operating system Administrator login for install and updates

NOV 22, 2016



Requirements for updating from 3.5 to 4.0

- Windows updates current is mandatory
- Microsoft .NET Framework 3.5 SP1
- Microsoft .NET Framework 4.0 Client Profile
- Windows operating system Administrator login for install and updates

Please **contact** Microquest for further information regarding network requirements/set-up.